



Howard Hughes Medical Institute
Research Laboratories

Stephen J. Elledge, Ph.D.
Investigator

Dear Colleagues,

Enclosed are several bacterial strains and a user manual for the Univector Plasmid Fusion System (UPS) that we have recently developed. The paper describing this system is currently being submitted. As you can see there are a number of recipient vectors included in this kit. Due to the large number of strains, we have aliquoted them into a sterile 96-well plate that is sealed. The strains have been stored frozen at -80°C in 9% DMSO prior to being sent to you. We have found that these strains can be successfully refrozen and stored after having been thawed for several days, so if you are not able to individually streak them out upon arrival you may be able to store them at -80°C and retrieve them later. You may also access the strains by poking a hole through the cover above the well of interest. A key for the identification of the strains by well position is included. Please remember that several strains bear no plasmids and must be grown with LB alone and the pUNI plasmid strains require kanamycin selection. The remaining strains are Amp resistant. Maps of the vectors are also included. The nomenclature for the vectors is pH for Host followed by a letter designating the organism the expression vector is designed for, B-bacterial, Y, yeast, I insect cells, M, mammalian. Each vector has a number to define the base vector and this is followed by the identity of the tag. For example, pHB1-HA3 is a bacterial vector that generates fusion proteins with a triple HA tag. pHB1-MYC3 has a triple MYC epitope tag. In addition there is a plasmid identification number for each plasmid that identifies that particular plasmid in our strain collection.

The plasmids we are sending are new and have been made by my lab and Dr. Wade Harper's lab. Dr. Harper has generously allowed us to distribute them prior to publication. We have shown that the UPS system works for bacterial, yeast, mammalian and insect cell expression systems as well as conventional vectors designed for the same purpose. While most of these plasmids have been tested, a few have not yet been used. We would be interested in hearing about any success stories or problems. Please address email to the univector mail box univ@bcm.tmc.edu. We are in the process of generating a number of vectors that will facilitate the conversion of your favorite vector into an epitope-tagged lox containing recipient vector. If you are interested in these, contact us in the future. Cre protein is available from commercial sources but we have found that making our own using the Gst-Cre fusion works much better in vitro.

If you happen to make new univector recipient vectors with novel properties, we would like to obtain a sample of the vector. If you could send your plasmid to us with a map after it has been successfully tested, we would appreciate that.

Sincerely,

Steve Elledge

Overview:

The Univector Plasmid-fusion System, UPS, is a new cloning method developed in our lab for rapid, efficient and systematic construction of recombinant DNA molecules. UPS employs the *cre-lox* site-specific recombination of bacteriophage P1 (Sternberg et al. 1981) to catalyze *in vitro* plasmid fusion events between the Univector, a plasmid containing the gene of interest, and host vectors containing the regulatory information. Plasmid fusion events are genetically selected and result in placement of the gene of interest under the control of novel regulatory elements. UPS eliminates the need for restriction enzymes, DNA ligase, and many *in vitro* manipulations required for subcloning. It greatly simplifies the process of designing and performing subclonings and saves researchers significant amounts of time and expense.

The Univector, as shown in Figure 1, contains a *loxP* site placed directly adjacent to the 5' end of the polylinker for insertion of the gene or cDNAs. The *loxP* site has a single open reading frame that is in frame with the ATG of the NdeI and NcoI sites of the polylinker. This facilitates the subsequent generation of protein fusions as noted below. Following the polylinker are bacterial and eukaryotic transcriptional terminators to facilitate 3' end formation of transcripts. Most importantly, the Univector has a conditional origin of replication derived from R6K γ that allows its propagation only in bacterial hosts expressing the *pir* gene. Finally, the Univector has the *neo* gene from *Tn5* for selection in bacteria.

Host vectors include the Col E1 origin of replication and the *bla* gene for propagation and selection in bacteria, a *loxP* site for plasmid fusion and a specific promoter residing upstream of and adjacent to the *loxP* site. In the case in which the pHOST plasmid also contains sequences for translational initiation, e.g. a GST tag, the pHOST *loxP* site is placed in the same reading frame as the *loxP* site on the Univector and the target gene must be cloned in the Univector in frame with the *loxP* open reading frame. Host vectors may also contain sequences responsible for propagation, selection and maintenance in organisms other than *E. coli*. Our lab is making a large collection of host vectors that are suitable for expression in bacteria, yeast, mammalian and insect cells. Some of these host vectors have been made and tested. They are listed in the table below (Figure 2). In general, host vectors are derivatives of existing vectors by insertion of an oligo linker containing a *loxP* or *loxH*¹ site in the polylinker of these vectors. For pHOST plasmids that also have epitope tags, the *loxP* site is inserted in frame with the coding sequence of the epitope tags.

Once the gene of interest is cloned in the Univector by conventional means, it can be rapidly transferred into multiple host vectors by UPS (Figure 3). The UPS reactions are carried out *in vitro* in the presence of purified GST-Cre fusion proteins. The GST-Cre fusion protein is easily produced and purified from *E. coli* in large quantities with glutathione-sepharose beads (see GST-Cre purification protocol). Plasmid fusions between the Univector and the pHOST plasmid via the *loxP* site create a plasmid dimer, thereby placing the target gene under the control of the regulatory sequences on the host vector. Selection for the recombinant products of UPS is achieved by selecting for kanamycin resistance after transformation into a *pir*-*E. coli* strain (standard cloning strains such as XL1 Blue, DH5 α , JM107, etc.) because the *neo* gene on the Univector can only be propagated when covalently linked to an origin of replication that is functional in a *pir*-background. UPS reactions can also be carried out *in vivo* but this method is slightly more complicated.

¹ *LoxP* consists of two 13 bp inverted repeats flanking an 8 bp central sequence. *loxH* is a mutant derivative of *loxP*, which has mutations at 3, 6, 9 nucleotide of the left repeat, creating three mismatches between the two repeats. The sequence is shown on the vector maps.

Cloning genes into the Univector

One of the main advantages of cloning by UPS is that, for each gene of interest, only one conventional subcloning is needed. Once the gene or cDNA is appropriately cloned into the Univector, it can be used to rapidly convert into any kind of recipient vectors to make expression constructs. Cloning by UPS is much simpler and faster than the conventional method. Not only does UPS avoid the use of restriction enzymes, DNA ligase and many *in vitro* manipulations for subcloning, it also saves researchers time in designing subcloning strategies. Below are the two key rules for cloning genes into the Univector.

1. The gene, or cDNA should be cloned in the same reading frame as the *loxP* site (Notice the *loxP* sequence contains only one open reading frame in the direction shown). Only then can it be used to seamlessly create expression constructs, transcriptional as well as translational fusions. There is no need to worry about reading frames or restriction site incompatibility. Everything should be automatically in frame because of the precise nature of the *cre-lox* site specific recombination.
2. It is recommended that one engineer a restriction site at the start codon, eg. NdeI or NcoI, such that only the ORF and no extra 5' untranslated sequence will be brought into the Univector. This will ease the construction of epitope-tagged constructs.

In vitro UPS reaction:

In general, the UPS reaction are carried out using 0.1-0.4 μg DNA each of the Univector and host vector in a volume of 20 μl in 1X Buffer S (50mM Tris.Cl pH 7.5, 10 mM MgCl_2 , 30 mM NaCl, 0.1 mg/ml BSA). The ratio of the two plasmids does not matter much but you will get more kanamycin resistant transformants if more Univector plasmid is used. We normally use 1 μl of the GST-Cre elution (0.1 to 0.2 μg) for each reaction. The reaction mixture is assembled on ice, incubated immediately at 37°C for 20 minutes and followed by 65°C for 5 minutes to inactivate the GST-Cre enzyme. However, the Cre inactivation is probably not necessary if you are going to immediately transform the mixture. Half of the reaction is used to transform heat-shock competent DH5 α or JM107 cells and transformants are plated onto LB/kanamycin (40 $\mu\text{g}/\text{ml}$) plates. Typically you will get 100 to 10,000 colonies. If you are using electroporation, 1 μl of the reaction added to 50 μl of cells is sufficient and will not arch. However, if one is performing a library conversion it is best to concentrate the DNA and remove the salt by precipitation or other methods prior to electroporation.

The efficiency of the GST-Cre mediated UPS reaction approaches 17%, as measured by the total number of Kn^{R} transformants vs the total number of Ap^{R} transformants, when equal molar amounts of the Univector and pHOST plasmid are used in the UPS reaction. By restriction analysis, 100% of the plasmids isolated from Kn^{R} transformants are correct fusion products via the *loxP* site. Most are dimers, with occasional trimers that in theory may result from two sequential fusion events or single fusion event between a pre-existing monomeric and a dimeric substrates. The production of trimers can be eliminated if gel purified monomeric supercoiled DNA is used suggesting the latter explanation is correct. Trimeric constructs usually express as well as dimers.

Occasionally we have found that UPS reaction give us exclusively trimers. In this case the Univector subclone was found to be a dimer. This occurs occasionally when too much DNA is used in ligation reactions (remember pUNI is puny). It is best to look at your Univector subclone uncut on a gel to determine if it is a monomer or dimer if it appears you are obtaining trimers from the UPS reaction.

RPM prepped DNA performs as well as or better than CsCl prepped DNA in our hands. We can routinely perform a dozen UPS reactions plus transformations in an hour. We normally examine only two Kn^{R} transformants from each plate and at least one of them will turn out to be a dimer. In most cases, both are dimers.

Restriction Analysis of the UPS products

To keep the digestion pattern analysis simple, it is preferable to use enzymes that cut only once or twice in each plasmid. Try to avoid enzymes that cut near the *loxP* site on both parental plasmids because the digestion pattern will not be informative. It is also recommended to digest both parental plasmids along with the UPS products to use as a reference. The digestion pattern for the dimeric UPS product should show certain bands missing and new bands appearing when compared to the parental plasmids. Restriction analysis of a trimer will produce one parental sized band in addition to the new bands.

Bacterial Strains (Relevant Genotype in Bold)


<u>Strain name</u>	<u>Genotype</u>
BW23473	Δ lac-169 rpoS(Am) robA1 creC510 hsdR514 endA recA1 uidA(Δ MluI):: pir⁺(wt) (This strain allows low copy pUNI replication, ~15 copies per cell)
BW23474	Δ lac-169 rpoS(Am) robA1 creC510 hsdR514 endA recA1 uidA(Δ MluI):: pir-116 (This strain allows high copy pUNI replication, ~250 copies per cell)
BUN10	hisG4 thr-1 leuB6 thi-1 lacY1 galK2 ara-14 xyl-5 mtl-1 argE3 rfbD1 mgl-51 kdgK51 Δ (gpt-proA)62 rpsL31 tsx33 supE44 recB21 recC22 sbcA23 hsdR::cat-pir-116(Cm^R) (This strain allows high copy pUNI replication and has the proper genetic background for homologous recombination and blue/white screening for lacO reconstitution. It is endA ⁺ so miniprep procedures using phenol extraction need to be used.)

GST-Cre purification protocol:

1. Grow a 500ml culture of the BL21 (pQL123 GST-Cre) strain to $OD_{600} = 0.8$ in LB/Amp (100 $\mu\text{g/ml}$) media.
2. Add IPTG to a final concentration of 0.4 mM. Induce for 3h at 37°C or 6h at room temperature. From our experience, GST-Cre prepared from the room temperature induction has a higher specific activity.
3. Pellet cells 5K for 10 min. You can freeze the cell pellet at -80°C at this point.
4. Resuspend the cell pellet in ice cold 10 ml NETN freshly supplemented with 1mM PMSF+ 5 $\mu\text{g/ml}$ leupeptin + 5 $\mu\text{g/ml}$ antipain and transfer to a clean 50 ml Falcon tube for sonication.
(NETN: 0.5% NP-40, 20mM TrisHCl pH 8.0, 100mM NaCl, 1mM EDTA)
5. Sonicate the 10 ml of cells at 4°C at level "2", 5 x 15 sec with an interval of 30 sec on ice between sonications. Try to avoid foaming which will denature the protein by keeping the sonicator tip well below the surface of the solution. The solution will become homogeneous and viscous during sonication.
6. Aliquot into six 2 ml eppendorf tubes, spin at 14,000 rpm for 15 min at 4°C to clear the cell extract. During this clearing step, prepare the glutathione-sepharose beads needed in step 7.
7. Prepare glutathione-sepharose beads (Pharmacia) in 200 μl per 2 ml eppendorf tubes. Wash three times with 1.5 ml NETN. Remove excess liquid, add the supernatant from step 6 to the beads and rotate at 4°C for one hour. The beads should be treated gently. We normally pellet them at 2500 rpm for 15 sec. using an eppendorf centrifuge.
8. Wash beads 4 times with 1.5 ml of ice cold NETN + protease inhibitors as in (4). Simply add the buffer, invert several times and spin at 2,500 rpm for 15 sec. This can be done at your bench for convenience.
9. Wash two times with wash buffer to replace the NETN prior to elution. Pellet the beads and remove all of the supernatant with thin-tip pipette tips (Phenix) that won't take up the beads.
(Wash Buffer: 100 mM Tris.Cl (pH 8.0), 120 mM NaCl, 1mM PMSF)
10. Elute three times by adding 150 μl elution buffer, rotating the tubes each time at 4°C for 45 minutes followed by a 15 sec spin at 2500 rpm. Collect and save supernatant each time. Typically you will get 0.1-0.2 μg GST-Cre protein per μl of elution buffer.
(Elution Buffer: 100 mM Tris.Cl (pH 8.0), 120 mM NaCl, 1mM PMSF, 20 mM Glutathione).
11. Aliquot the GST-Cre elution into 80 μl aliquots or smaller and store at -80°C . Check the quantity and quality of the GST-Cre prep by running 5 μl of the elution on a 10% SDS PAGE gel. Once a tube is taken out, it can be stored at 4°C for up to a month. However, the enzymatic activity of GST-Cre will decrease with time.

Figure 1

A



Epitope Tag Sequences

1. HA

```
Nco I                               Nhe I                               Nde I
|                                   |                                   |
CC ATG GCT TAC CCA TAC GAT GTT CCA GAT TAC GCT AGC TTG GGT GGT CAT ATG
M  A  Y  P  Y  D  V  P  D  Y  A  S  L  G  G  H  M
3/1                                33/11
```

2. HA3

```
Nco I
|
CC ATG GCA GGT TAC CCA TAC GAC GTT CCT GAC TAT GCG TCA CTC TAC CCC TAT GAC GTA CCG
M  A  G  Y  P  Y  D  V  P  D  Y  A  S  L  Y  P  Y  D  V  P
6/1                                36/11

Nsi I                               Xba I
|                                   |
GAT TAT GCA TCC CTA TAT CCG TAT GAT GTT CCA GAT TAC GCT TCT CTA CGT TCC TCT AGA
D  Y  A  S  L  Y  P  Y  D  V  P  D  Y  A  S  L  R  S  S  R
66/21                                96/31

Nde I
|
GGC GTC CAC CAT ATG
G  V  H  H  M
126/41
```

3. MYC3

```
NcoI                               Eag I                               NotI
|                                   ||                               ||
CC ATG GGA TTC GAG CTA TGC GGC CGC TCT GAG CAA AAG CTC ATT TCT GAA GAG GAC TTG AAT
M  G  F  E  L  C  G  R  S  E  Q  K  L  I  S  E  E  D  L  N
3/1                                33/11

Bcl I
|
GGA GAA CAG AAA TTG ATC AGT GAG GAA GAC CTC AAC GGT GAG CAG AAG TTA ATA TCC GAG
G  E  Q  K  L  I  S  E  E  D  L  N  G  E  Q  K  L  I  S  E
63/21                                93/31

Eag I                               Not I                               Spe I                               Apa I                               Nde I
||                               ||                               |                               |                               |
GAG GAT CTT AAT AGT AGC GGC CGC ACT AGT CTC GAA GGG GGG CCC CAT ATG
E  D  L  N  S  S  G  R  T  S  L  E  G  G  P  H  M
123/41                                153/51
```

4. MYC9

```
Nco I
|
CC ATG GCT AGT (GGT GAA CAA AAG TTG ATT TCT GAA GAA GAT TTG AAC GGT AGG CAA AAG CTA
M  A  S  G  E  Q  K  L  I  S  E  E  D  L  N  G  R  Q  K  L
3/1                                33/11

BamH I
|
ATC TCC GAG GAA GAC TTC AAC GGT GAA CAA AAA TTA ATC TCA GAA GAA GAC TTG AAC GGA
I  S  E  E  D  F  N  G  E  Q  K  L  I  S  E  E  D  L  N  G
63/21                                93/31

Xba I                               Spe I
|                                   |
TCC TCT AGA) 3 ACT AGT
S  S  R  T  S
123/41
```

Note: The sequence in the parenthesis repeats three times.

Epitope Tag Sequences

1. HA

```
Nco I                               Nhe I                               Nde I
|                                   |                                   |
CC ATG GCT TAC CCA TAC GAT GTT CCA GAT TAC GCT AGC TTG GGT GGT CAT ATG
M  A  Y  P  Y  D  V  P  D  Y  A  S  L  G  G  H  M
3/1                                33/11
```

2. HA3

```
Nco I
|
CC ATG GCA GGT TAC CCA TAC GAC GTT CCT GAC TAT GCG TCA CTC TAC CCC TAT GAC GTA CCG
M  A  G  Y  P  Y  D  V  P  D  Y  A  S  L  Y  P  Y  D  V  P
6/1                                36/11

Nsi I                               Xba I
|                                   |
GAT TAT GCA TCC CTA TAT CCG TAT GAT GTT CCA GAT TAC GCT TCT CTA CGT TCC TCT AGA
D  Y  A  S  L  Y  P  Y  D  V  P  D  Y  A  S  L  R  S  S  R
66/21                                96/31

Nde I
|
GGC GTC CAC CAT ATG
G  V  H  H  M
126/41
```

3. MYC3

```
NcoI                               Eag I                               NotI
|                                   ||                               ||
CC ATG GGA TTC GAG CTA TGC GGC CGC TCT GAG CAA AAG CTC ATT TCT GAA GAG GAC TTG AAT
M  G  F  E  L  C  G  R  S  E  Q  K  L  I  S  E  E  D  L  N
3/1                                33/11

Bcl I
|
GGA GAA CAG AAA TTG ATC AGT GAG GAA GAC CTC AAC GGT GAG CAG AAG TTA ATA TCC GAG
G  E  Q  K  L  I  S  E  E  D  L  N  G  E  Q  K  L  I  S  E
63/21                                93/31

Eag I                               Not I                               Spe I                               Apa I                               Nde I
||                               ||                               |                               |                               |
GAG GAT CTT AAT AGT AGC GGC CGC ACT AGT CTC GAA GGG GGG CCC CAT ATG
E  D  L  N  S  S  G  R  T  S  L  E  G  G  P  H  M
123/41                                153/51
```

4. MYC9

```
Nco I
|
CC ATG GCT AGT (GGT GAA CAA AAG TTG ATT TCT GAA GAA GAT TTG AAC GGT AGG CAA AAG CTA
M  A  S  G  E  Q  K  L  I  S  E  E  E  D  L  N  G  R  Q  K  L
3/1                                33/11

BamH I
|
ATC TCC GAG GAA GAC TTC AAC GGT GAA CAA AAA TTA ATC TCA GAA GAA GAC TTG AAC GGA
I  S  E  E  D  F  N  G  E  Q  K  L  I  S  E  E  D  L  N  G
63/21                                93/31

Xba I                               Spe I
|                                   |
TCC TCT AGA) 3 ACT AGT
S  S  R  T  S
123/41
```

Note: The sequence in the parenthesis repeats three times.

5. GST

Xba I

|
TCTAG ATG CTA GGT TAT TGG AAA ATT AAG GGC CTT GTG CAA CCC ACT CGA CTT CTT TTG GAA TAT
M L G Y W K I K G L V Q P T R L L L E Y
6/1 36/11

CTT GAA GAA AAA TAT GAA GAG CAT TTG TAT GAG CGC GAT GAA GGT GAT AAA TGG CGA AAC
L E E K Y E E H L Y E R D E G D K W R N
66/21 96/31

AAA AAG TTT GAA TTG GGT TTG GAG TTT CCC AAT CTT CCT TAT TAT ATT GAT GGT GAT GTT
K K F E L G L E F P N L P Y Y I D G D V
126/41 156/51

AAA TTA ACA CAG TCT ATG GCC ATC ATA CGT TAT ATA GCT GAC AAG CAC AAC ATG TTG GGT
K L T Q S M A I I R Y I A D K H N M L G
186/61 216/71

GGT TGT CCA AAA GAG CGT GCA GAG ATT TCA ATG CTT GAA GGA GCG GTT TTG GAT ATT AGA
G C P K E R A E I S M L E G A V L D I R
246/81 276/91

TAC GGT GTT TCG AGA ATT GCA TAT AGT AAA GAC TTT GAA ACT CTC AAA GTT GAT TTT CTT
Y G V S R I A Y S K D F E T L K V D F L
306/101 336/111

AGC AAG CTA CCT GAA ATG CTG AAA ATG TTC GAA GAT CGT TTA TGT CAT AAA ACA TAT TTA
S K L P E M L K M F E D R L C H K T Y L
366/121 396/131

AAT GGT GAT CAT GTA ACC CAT CCT GAC TTC ATG TTG TAT GAC GCT CTT GAT GTT GTT TTA
N G D H V T H P D F M L Y D A L D V V L
426/141 456/151

TAC ATG GAC CCA ATG TGC CTG GAT GCG TTC CCA AAA TTA GTT TGT TTT AAA AAA CGT ATT
Y M D P M C L D A F P K L V C F K K R I
486/161 516/171

Sca I

|
GAA GCT ATC CCA CAA ATT GAT AAG TAC TTG AAA TCC AGC AAG TAT ATA GCA TGG CCT TTG
E A I P Q I D K Y L K S S K Y I A W P L
546/181 576/191

CAG GGC TGG CAA GCC ACG TTT GGT GGT GGC GAC CAT CCT CCA AAA TCG GAT CTG GTT CCG
Q G W Q A T F G G G D H P P K S D L V P
606/201 636/211

Nde I Nco I

| |
CGT GGA TCT CGT CGT GCA TCT GTT GGA TCG CAT ATG CCC ATG G
R G S R R A S V G S H M P M
726/241 756/251

6. His6

```
Nco I                                     Nde I
|                                         |
CC ATG GGC AGC AGC CAT CAT CAT CAT CAT CAC AGC AGC GGC CTG GTG CCG CGC GGC AGC CAT ATG
M   G   S   S   H   H   H   H   H   H   S   S   G   L   V   P   R   G   S   H   M
3/1                                     33/11                                     63/21
```

7. FLAG


```
Nco I                                     Nde I
|                                         |
CC ATG GAC TAC AAA GAC GAT GAC GAC AAG CAT ATG
M   D   Y   K   D   D   D   D   K   H   M
3/1                                     33/11
```

8. GFP

9. EE

```
Nco I                                     Aat II  Nde I
|                                         |         |
CC ATG GAA GAA GAG GAA TAC ATG CCA ATG GAG GAC GTC CAT ATG
M   E   E   E   E   Y   M   P   M   E   D   V   H   M
3/1                                     33/11
```

pTag-lox



Bluescript II KS backbone

Other tags: Flag, HA, HA3, His6, Myc3, Myc9, GST

The pTAG vectors are designed to aid in the construction of pHOST vectors containing commonly used epitope and affinity tags. Each cassette contains a tag followed by a kanamycin resistance gene, neo from Tn5, flanked by lox sites. This cassette can be excised using restriction sites outside of the cassette such as SmaI and ligated into the vector of choice by selecting for kanamycin resistance and by identifying clones containing the proper orientation. The neo gene can then be removed by cre-lox recombination in vitro using the Gst-Cre enzyme, or in vivo by transformation into a cre expressing strain such as BNN132. If you use BNN132 it is best to streak out transformants once or retransform the plasmid into a second strain to ensure there are no non-recombinant plasmids contaminating the prep because recombination in BNN132 is not 100% efficient. Plasmids that have undergone recombination between lox sites can now be used as recipient expression vectors for UPS.

pHOST (pH) Vectors for UPS

Bacterial pH vectors (pHB1-99):

Plasmid Identification Numbers

Vectors	Origin	loxH	HA	HA3	MYC3	MYC9	GST	His6	FLAG	GFP	EE
pHB1	pETX (T7)		p1201	p1202	p1203	p1204					
pHB2	pGEX-2T (tac)						p1205				
pHB3	pET15b (T7)							p1206	p1207		

Insect cell pH vectors (pHI100-199):

Plasmid Identification Numbers

Vectors	Origin	loxP	HA	HA3	MYC3	MYC9	GST	His6	FLAG	GFP	EE
pHI100	pVL1392	p1208		p1210	p1211		p1212	p1213	p1214		p1215

Mammalian pH vectors (pHM200-299):

Plasmid Identification Numbers

Vectors	Origin	loxH	HA	HA3	MYC3	MYC9	GST	His6	FLAG	GFP	EE
pHM200	pcDNA3.1			p1225	p1226						
pHM201	pFLAG-CMV2								p1227		

Yeast (*S. cerevisiae*) pH vectors (pHY300-):

Plasmid Identification Numbers

Vectors	Origin	Features	loxH	HA	HA3	MYC3	MYC9	GST	His6	FLAG	GFP	EE
pHY313	pRS413	CEN HIS GAL										
pHY314	pRS414	CEN TRP GAL	p1216			p1217	p1218	p1219				
pHY315	pRS415	CEN LEU GAL										
pHY316	pRS416	CEN URA GAL	p1220									
pHY323	pRS423	2 μ HIS GAL										
pHY324	pRS424	2 μ TRP GAL	p1221									
pHY325	pRS425	2 μ LEU GAL	p1222									
pHY326	pRS426	2 μ URA GAL	p1223					p1224				

Plate position	Vector Names	Plasmid ID #	Base Vector	Description
2A				
2B	pHB1-HA	p1201	PETX (T7)	<i>E. coli</i> expression vector for the generation of HA-tagged fusion proteins
2C	pHB1-HA3	p1202	PETX (T7)	<i>E. coli</i> expression vector for the generation of 3xHA-tagged fusion proteins
2D	pHB1-MYC3	p1203	PETX (T7)	<i>E. coli</i> expression vector for the generation of 3xMYC-tagged fusion proteins
2E	pHB1-MYC9	p1204	PETX (T7)	<i>E. coli</i> expression vector for the generation of 9xMYC-tagged fusion proteins
2F	pHB2-GST	p1205	pGEX-2T (tac)	<i>E. coli</i> expression vector for the generation of GST-tagged fusion proteins
2G	pHB3-His6	p1206	PET15b (T7)	<i>E. coli</i> expression vector for the generation of His6-tagged fusion proteins
2H				
3A				
3B	pHB3-FLAG	p1207	PET15b (T7)	<i>E. coli</i> expression vector for the generation of FLAG-tagged fusion proteins
3C	pHI100-loxP	p1208	pVL1392	Baculovirus expression vector for the expression of untagged proteins
3D				
3E	pHI100-HA3	p1210	pVL1392	Baculovirus expression vector for the expression of 3xHA-tagged fusion proteins
3F	pHI100-MYC3	p1211	pVL1392	Baculovirus expression vector for the expression of 3xMYC-tagged fusion proteins
3G	pHI100-GST	p1212	pVL1392	Baculovirus expression vector for the expression of GST-tagged fusion proteins
3H				
4A				
4B	pHI100-His6	p1213	pVL1392	Baculovirus expression vector for the expression of His6-tagged fusion proteins
4C	pHI100-FLAG	p1214	pVL1392	Baculovirus expression vector for the expression of FLAG-tagged fusion proteins
4D	pHI100-EE	p1215	pVL1392	Baculovirus expression vector for the expression of GluGlu-tagged fusion proteins
4E	pHY314-loxH	p1216	pRSS414	<i>S. cerevisiae</i> CEN TRP1 GAL expression vector for the expression of untagged proteins

Amp
Amp
Amp
Amp
Amp
Amp

Amp
Amp
Amp

Amp
Amp

Amp
Amp

Amp
Amp

Plate position	Vector Names	Plasmid ID #	Base Vector	Description
4F	pHY314-HA3	p1217	pRS414	<i>S. cerevisiae</i> CEN TRP1 GAL expression vector for the expression of 3xHA-tagged fusion proteins
4G	pHY314-MYC3	p1218	pRS414	<i>S. cerevisiae</i> CEN TRP1 GAL expression vector for the expression of 3xMYC-tagged fusion proteins
4H				
5A				
5B	pHY314-MYC9	p1219	pRS414	<i>S. cerevisiae</i> CEN TRP1 GAL expression vector for the expression of 9xMYC-tagged fusion proteins
5C	pHY316-loxH	p1220	pRS416	<i>S. cerevisiae</i> CEN URA3 GAL expression vector for the expression of untagged proteins
5D	pHY324-loxH	p1221	pRS424	<i>S. cerevisiae</i> 2 μ TRP GAL expression vector for the expression of untagged proteins
5E	pHY325-loxH	p1222	pRS425	<i>S. cerevisiae</i> 2 μ LEU GAL expression vector for the expression of untagged proteins
5F	pHY326-loxH	p1223	pRS426	<i>S. cerevisiae</i> 2 μ URA GAL expression vector for the expression of untagged proteins
5G	pHY326-MYC9	p1224	pRS426	<i>S. cerevisiae</i> 2 μ URA GAL expression vector for the expression of untagged proteins
5H				
6A				
6B	pHM200-HA3	p1225	pcDNA3.1	Mammalian expression vector for the expression of 3xHA tagged fusion proteins
6C	pHM200-MYC3	p1226	pcDNA3.1	Mammalian expression vector for the expression of 3xMYC tagged fusion proteins
6D	pHM201-FLAG	p1227	pFLAG-CMV2	Mammalian expression vector for the expression of FLAG tagged fusion proteins
6E	pUNI 10	p1228		Original Univector
6F	pUNI 15	p1229		A Univector derivative with an fl ori for ss DNA production

Kan Kan

Amp Amp

Amp Amp

Amp Amp

Amp Amp

Amp Amp

Amp Amp

Amp Amp

Amp Amp

Amp Amp


Amp Amp

Amp Amp

Plate position	Vector Names	Plasmid ID #	Base Vector	Description
6G	pUNI 20	p1230		A Res site and a Rs site were inserted into pUNI 15.
6H				
7A				
7B	PACT2lox	p1231		A 2 μ LEU2 ADH-GAL4 (AD) fusion vector for the 2 hybrid system
7C	PAS2lox	p1232		A 2 μ TRP1 ADH-GAL4 (DBD) fusion vector for the 2 hybrid system
7D	pDAB1lox	p1233		A CEN TRP1 ADH-GAL4 (DBD) fusion vector for the 2 hybrid system
7E	PSOC4lox	p1234		A 2 μ URA3 ADH-hSOS vector for SOS fusion for the SRS system
7F	pMS-TRPlox	p1235		A 2 μ TRP1 GAL-Myristoylation vector for the SRS system
7G				
7H				
8A				
8B	BW23473			Bacterial strain (hsdR endA1 recA1 pir+) for low copy propagation of pUNI clones
8C	BW23474			Bacterial strain (hsdR endA1 recA1 pir-116) for high copy propagation of pUNI clones
8D	BUN10			Bacterial strain (hsdR ^{RM} ::CAT recBC sbcA pir-116) for homologous recombination of pUNI clones
8E	pGST-cre	p1236		Bacterial strain containing a plasmid (pGST-cre) overproducing GST-cre for in vitro cre-lox reactions.
8F	pTag-FLAGlox	p1237		pTag vector for transferring an epitope tag (FLAG) and lox into a host vector.
8G	pTag-HAlox	p1238		pTag vector for transferring an epitope tag (HA) and lox into a host vector.
8H				
9A				
9B	pTag-HA3lox	p1239		pTag vector for transferring an epitope tag (HA3) and lox into a host vector.
9C	pTag-MYC3lox	p1240		pTag vector for transferring an epitope tag (MYC3) and lox into a host vector.

Kan
Amp
Amp
Amp
Amp
Amp
Amp

pHB1-HA /HA3 /MYC3 /MYC9




pHB1 polylinker sequence: (only unique sites are shown)

Nco I Nde I LoxP
 CC ATG G (epitope tag sequence) CAT ATG ATA ACT TCG TAT AGC ATA

BamH I
 CAT TAT ACG AAG TTA TCG ATA GGA TCC


pHB2-GST:



pHB2-GST polylinker sequence (unique sites are shown)

Nde I Nco I loxP
 CAT ATG CCC ATG GCT ATA ACT TCG TAT AGC ATA CAT TAT ACG AAG

BamH I EcoR I
 TTA TGG ATC CGA ATT C




pHB3-His6/FLAG polylinker sequence:
 (only unique sites are shown)

Nco I Nde I LoxP
 CC ATG G (epitope tag sequence) CAT ATG ATA ACT TCG TAT AGC ATA

BamH I
 CAT TAT ACG AAG TTA TCG ATA GGA TCC

pHM200-HA3/ MYC3
(pcDNA 3.1 delta NEO)




pHM200-HA3/ MYC3 polylinker sequence:
(unique sites are shown in bold)

Nco I Nde I LoxP
CC ATG G (epitope tag sequence) CAT ATG ATA ACT TCG TAT AGC ATA CAT

BamH I
TAT ACG AAG TTA TCG ATA GGA TCC

pHI201-FLAG
(pFLAG-CMV-2)

pHM201-FLAG





pHI201-FLAG polylinker sequence: (only unique sites are shown)

Hind III Not I Aat II loxP
AAG CTT GCG GCC GCG GAC GTC ATA ACT TCG TAT AGC ATA CAT TAT

BamH I
ACG AAG TTA TAC GAA GTT ATG GAT CCC GGG
Sma I

More information about pFLAG-CMV-2 may be obtained from Sigma, 1-888-397-4542



pHY314-loxH polylinker sequence:
(unique sites are shown in bold)

LoxH

GAA TTC GAC GTC ATT ACC TCA TAT AGC ATA CAT TAT ACG AAG TTA

* * *

Sal I Pst I Sma I

TTG TCG ACG GTA TCG ATA AGC TTG ATA TCG AAT TCC TGC AGC CCG

BamH I Spe I Not I Sac II Sac I

GGG GAT CCA CTA GTT CTA GAG **CGG CCG** CCA CCG CGG TGG AGC TC

Eag I

* marks the positions of the nucleotide substitutions in loxH

pHY314-MYC9 polylinker sequence:
(unique sites are shown in bold)

LoxP

CC ATG G (MYC9 tag sequence) AGT ATG ATA ACT TCG TAT AGC ATA

* * *

Nco I Cla I BamH I Spe I Xba I Not I


CAT TAT ACG AAG TTA TCG ATA GGA TCC ACT AGT TCT AGA **GCG GCC**

Eag I


Sac II Sac I

GCC ACC GCG GTG GAG CTC

pHY314-HA3/ MYC3
(CEN TRP GAL)



pHY316-loxH
(CEN URA GAL)



pHY314-HA3/ MYC3 polylinker sequence:
(unique sites are shown in bold)

LoxH

CC ATG G (epitope tag sequence) CAT ATG ATA ACT TCG TAT AGC ATA

* * *

Nco I Nde I LoxP

CC ATG G (epitope tag sequence) CAT ATG ATA ACT TCG TAT AGC ATA

Sal I Pst I Sma I

TTG TCG ACG GTA TCG ATA AGC TTG ATA TCG AAT TCC TGC AGC CCG

BamH I Spe I Not I

GGG GAT CCA CTA GTT CTA GAG **CGG CCG** CCA CCG CGG TGG AGC TC

Eag I

Sac II Sac I

GCC ACC GCG GTG GAG CTC

pHY316-loxH polylinker sequence:
(unique sites are shown in bold)

LoxH

GAA TTC GAC GTC ATT ACC TCA TAT AGC ATA CAT TAT ACG AAG TTA

* * *

Sal I Pst I Sma I


TTG TCG ACG GTA TCG ATA AGC TTG ATA TCG AAT TCC TGC AGC CCG

BamH I Spe I Not I

GGG GAT CCA CTA GTT CTA GAG **CGG CCG** CCA CCG CGG TGG AGC TC

Eag I

* marks the positions of the nucleotide substitutions in loxH




pHY326-loxH polylinker sequence:
(unique sites are shown in bold)

Aat II LoxH
GAA TTC GAC GTC ATT ACC TCA TAT AGC ATA CAT TAT ACG AAG TTA
* * *
Sal I Hind III Pst I Sma I
TTG TCG ACG GTA TCG ATA AGC TTG ATA TCG AAT TCC TGC AGC CCG
BamH I Spe I Not I Sac II Sac I
GGG GAT CCA CTA GTT CTA GAG **CGG CCG** CCA CCG CGG TGG AGC TC
Eag I


* marks the positions of the nucleotide substitutions in loxH

pHY326-MYC9
(2μ URA GAL)



pHY326-MYC9 polylinker sequence:
(unique sites are shown in bold)

Nco I LoxP
CC ATG G (MYC9 tag sequence) AGT ATG ATA ACT TCG TAT AGC ATA CAT
Cla I Hind III
TAT ACG AAG TTA TCG ATA AGC TTG ATA TCG AAT TCC TGC AGC CCG
Spe I Not I Sac II Sac I
TGG GAT CCA CTA GTT CTA GAG **CGG CCG** CCA CCG CGG TGG AG CTC
Eag I




pHY326-loxH polylinker sequence:
(unique sites are shown in bold)

Aat II LoxH
GAA TTC GAC GTC ATT ACC TCA TAT AGC ATA CAT TAT ACG AAG TTA
* * *
Sal I Cla I Hind III Sma I
TTG TCG ACG GTA TCG ATA AGC TTG ATA TCG AAT TCC TGC AGC CCG
BamH I Spe I Not I Sac II Sac I
GGG GAT CCA CTA GTT CTA GAG **CGG CCG** CCA CCG CGG TGG AGC TC
Eag I

* marks the positions of the nucleotide substitutions in loxH


pHY324-loxH
(2μ TRP GAL)



pHY324-loxH polylinker sequence:
(unique sites are shown in bold)

Aat II LoxH
GAA TTC GAC GTC ATT ACC TCA TAT AGC ATACAT TAT ACG AAG
Sal I Cla I Hind III
TTA TTG TCG ACG GTA TCG ATA AGC TTG ATA TCG AAT TCC TGC
Sma I BamH I Spe I Not I Sac II
AGC CCG GGG GAT CCA CTA GTT CTA GAG **CGG CCG** CCA CCG CGG
Sac I Eag I
TGG AGC TC


* marks the positions of the nucleotide substitutions in loxH



pACT2lox polylinker

LOX


CAT ATG GCC ATG GTA ACT TCG TAT AGC ATA CAT TAT ACG AAG TTA
 NdeI NcoI
 GTC GAG AGA TCT
 Bgl II



pAS2lox polylinker

LOX



CAT ATG GCC ATG GTA ACT TCG TAT AGC ATA CAT TAT ACG AAG TTA
 NdeI NcoI
 GTC GAC
 Sal I



pDAB1lox polylinker

LOX

CAT ATG GCC ATG GTA ACT TCG TAT AGC ATA CAT TAT ACG
 NdeI NcoI
 AAG TTA GTC GAC
 Sal I




pUNI polylinker sequence. Unique restriction sites are in bold.

(401) NotI KpnI LOX
 GC GGC CGC GGT ACC ATA ACT TCG TAT AGC ATA CAT TAT ACG AAG TTA TCT
EcoRI SmaI XhoI NdeI NcoI BamHI NotI
 GGA ATT CCC CGG GCT CGA GAA CAT ATG GCC ATG GGG ATC CGC GGC CGC
HpaI SalI SacI
 AAT TGT TAA CAG ATC CGT CGA CGA GCT CGC TA (530)

pUNI 15 Polylinker Sequence

(37) KpnI LOX EcoRI
 GGT ACC ATA ACT TCG TAT AGC ATA CAT TAT ACG AAG TTA TCT GGA ATT
SmaI XhoI NdeI NcoI BamHI NotI MunI
 CCC CGG GCT CGA GAA CAT ATG GCC ATG GGG ATC CGC GGC CGC AAT
HpaI SalI SacI (156)
 TGT TAA CAG ATC CGT CGA CGA GCT CGC



pUNI 20 Polylinker Sequence

(157) KpnI LOX EcoRI
 GGT ACC ATA ACT TCG TAT AGC ATA CAT TAT ACG AAG TTA TCT GGA ATT
SmaI XhoI NdeI NcoI BamHI NotI
 CCC CGG GCT CGA GAA CAT ATG GCC ATG GGG ATC CGC GGC CGC (246)